

ASPEO®



BONE COLLECTOR

EN - INSTRUCTIONS FOR USE

Anthogyr
PRIME MOVER IN IMPLANTOLOGY

TABLE OF CONTENTS

- I. MEANING OF SYMBOLS USED
- II. GENERAL INFORMATION
- III. DESCRIPTION
- IV. RECOMMENDATIONS
- V. IMPORTANT
- VI. ASSEMBLING / DISMANTLING
- VII. SURGERY PROCEDURE
- VIII. HYGIENE AND MAINTENANCE
- IX. REPAIR
- X. ACCESSORIES
- XI. WARRANTIES

I. MEANING OF SYMBOLS USED



These **MDs** comply with the Community Directive 93/42/CEE as amended by the 2007/47/CEE directive.



Caution: it is necessary to refer to the precautions of use for any important information related to safety



Do not reuse



Sterilised using irradiation



Not sterile



Wear gloves



Sterilisation by autoclave, at the specified temperature



Use by date



Catalog reference number



Batch code



Serial number



Do not use if package is damaged



Keep dry



Keep away from sunlight



Manufacturer




Date of manufacture

II. GENERAL INFORMATION

Aspeo® is to be used by practitioners who have completed the relevant practical and theoretical trainings. Through a systematic use of the Aspeo® bone collector, surgeries can be made more secure and optimised. Endogenous bone can be made available for every surgery to make up for the lack of bone or to better deal with the aesthetic aspect.

III. DESCRIPTION

The system is constituted by a main body with a detachable front section, a removable front canula, a cross tube to be plugged to the chair suction system, a movable piston, all parts in stainless steel, to be completed with an inner cylindrical filter, made of polyethylene high density, restrained at one end by the piston.



Each package contains:

- a main body (Stainless steel)
- a front section (Stainless steel)
- a front canula (Stainless steel)
- a piston fitted with a removable button. (Stainless steel)
- an adaptor (Plastic)
- three single use sterile filters (Plastic)

The device is to be plugged to the suction systems by means of a connector of inner Ø6,5mm or by using an adaptor (Ref.12041) for Ø11mm and Ø16mm. The Aspeo® system comes with 3 disposable filters. It must be used **EXCLUSIVELY** with one of the filters of batches Ref. 12006 or 12012 supplied by Anthogyr. Each filter is individually packaged in a sealed pouch and Gamma Ray sterilised. Only unaltered packaging can guarantee that the filter's sterility is maintained.

⚠ Used filters or filters featuring defective packaging should be immediately discarded. **ONLY** Anthogyr's filters Ref. 12006 or 12012 should be used with the Aspeo® system.

IV. RECOMMENDATIONS

The Aspeo® system is designed to gather bone fragments which can be used to perform autogenous bone fillings may it be deemed necessary during the implant surgery. It is particularly recommended for implant and periodontal surgery. Bone fragments harvested and mixed with the patient's blood may, for example, be used as space maintainers or to fill in a bone defect in order to facilitate rehabilitation during reconstruction procedures.

V. IMPORTANT

- ⚠** • The Aspeo® system must be cleaned, disinfected and sterilised before its first use and after each utilisation.
- ⚠** • Proceed stringently when working with the Aspeo® at the sampling

site and always use a second suction system for the saliva.

- Care should be taken to ensure that the harvested bone fragments are not dried by the continuous airflow going through the filter; when the device is not in use, shut off the suction airflow.

- Handle the harvested bone fragments with great care, complying with well established surgical protocol.

- Aspeo® may be used several times during a same intervention to harvest a significant volume of bone from/on the same patient.



- Connect the device to a CE medical suction system compliant with the EN ISO 11607 standard.



- Maintain copious irrigation to avoid burning the bone.

- If the bone is not placed immediately, it is recommended to deposit it in a sterile cup containing a sterile saline solution until ready to put in place.

- It is recommended to carry out a preoperative rinsing of the mouth.

- It is recommended to use a mouthwash for several days after the procedure to avoid bacterial growth in the mouth.



- Prior to use the Aspeo®, after its complete mounting and plugging onto the chair's suction system, make sure to pull the piston back until it clicks into its rear position. Take care not to engage the piston when the front canula is in position. This could result in ejection of the canula and of the bone fragments.

- In the event of a significant drop in suction, clean all the tube sections: the front canula, the filter and tube.



- All parts should be taken apart and cleaned immediately after each intervention.



- The filter is DISPOSABLE. It is to be discarded after use and eliminated through the waste disposal network.

Single-use devices : do not reuse and / or re-sterilise (potential risks of contamination and modification of functional surfaces).


VI. ASSEMBLING / DISMANTLING

To dismantle :

- 1 Disconnect Aspeo® from the suction line.

- 2 Disconnect the adaptor

- 3 Remove the front canula of the bone collector.





- 4 Unscrew the front section.

- 5 Remove the filter.

- 6 Unscrew the piston button.



- 7 Remove the piston by pushing it through the Aspeo® body.



Insert:

- 1 Insert the piston (without its button) from left to right into the body and pull it out entirely.

- 2 Put the filter in position in its housing.





- 3 Screw the front section on fully.

- 4 Screw on the piston button.

- 5 Pull the piston back until you hear a locking "click".

- 6 Replace the front canula.




VII. SURGERY PROCEDURE

Before use, ensure that the product is assembled and the piston is working correctly. The front canula should be set into its housing. The piston should be in the rear position.

- Connect the Aspeo® to the suction system.

- Keep the front canula close to the implant site while drilling.




TO DISPOSE OF THE BONE FRAGMENTS:

- Remove the front canula.

- Engage the piston by pushing it forward as on a syringe.

Collected bone fragments may then be deposited in a cup or directly placed on the site to be filled.




TO RESUME BONE FRAGMENT HARVESTING:

HARVESTING:

- Pull the piston back to its initial rear position. A "click" will confirm its correct positioning.

- Replace the front canula.

- Reconnect the Aspeo® to the suction system.



VIII. HYGIENE AND MAINTENANCE

The re-sterilisation of the reusable medical devices must be done by a properly trained and protected staff, while respecting regulations in place. The re-sterilisation protocol must be appropriate to match with the risks of infection.

For each device used please refer to the manufacturer's instructions. More particularly, please respect the concentration levels, the duration of the exposures, the liquid solutions replacement and the products' life-time. Never mix the products. Please dispose of the products as per instructed in the notices.

⚠ Use only products intended for the maintenance of the medico-surgical material.

Proscribe the antiseptics which are intended to be used for the skin and the mucous membranes.

⚠ Only use products compatible with the stainless steel.

① Appropriate protective clothes and gloves must be used in order to prevent any risk of infection and injury.

PRE-DISINFECTION

Pre-disinfection is to be conducted immediately after the operation.

By soaking:

- Use only detergent or disinfectant solutions with a neutral or slightly alkaline pH. It is not advised to use products which may fix X proteins to the unit (alcohol, aldehydes...).
- Carefully rinse and dry every part.

Or : thermo-disinfection (washer disinfectant), immediately after the surgery.

- Place carefully each part.
- The thermo-disinfection cycle lasts 10 minutes at 93°C (203°F). Follow the manufacturer's instructions for use.
- Use only low alkaline cleaning products with pH is between 2.5 and 9.
- Respect imperatively the drying cycle.

⚠ Drain the possible water residues with pressured medical quality air : Wear glasses and a mask.

CLEANING

- Use only detergent or disinfectant solutions with a neutral or slightly alkaline pH. It is not advised to use products which may fix X proteins to the unit (eg. Alcohol, aldehydes...). Immerse the device completely while respecting the duration of the cycle.
- This dismounted device can be immersed 60 minutes in sodium hypochlorite 1M or 2M. Do not use an aluminium container.

⚠ The use of sodium hypochlorite (bleach) is to be proscribed: Important risk of corrosion.

Use a low frequency ultrasonic tank (25 to 50 KHz) by heating the solution to the maximum of 45° or brush with a soft brush and clean each element very carefully (a single use brush is recommended).

⚠ Contacts of parts against each other or with the tank can involve strains.

⚠ Make sure, for the use of an ultrasonic tank, that the solution of detergent disinfecting is compatible with this process.

Rinse the parts with softened water, then with a demineralized water in a low frequency ultrasonic tank (25 to 50KHz) for 10 minutes at 60°C.

Dry each pieces carefully, with soft and non-fluffy sterile pads, with the possible use of pressured medical air.

Check if the cleaning succeeded in eliminating any visible stain, and if not, redo the cleaning procedure.

STERILISATION

The device must be sterilised prior and after the first and every use.

The device and its accessories can be sterilised at 135°C (275°F) at a pressure of 2.13 bars **only** in a fractionated vacuum steam autoclave for a minimum of 20 minutes (sterilisation maintenance period). **The single use filter is delivered sterile.**

We strongly advise to use class B autoclaves.

Any other sterilisation method is to be proscribed.

ⓘ Read the autoclave manufacturer's instructions of use. Respect a certain distance between each bag and do not overload the autoclave.

- Before sterilisation, the instruments must be pre-disinfected, cleaned, and tested.
- Reassemble all the parts of the Aspeo® but for the front cannula, without any filters.
- Check that the device does not have any point of corrosion or cracks and control its correct functioning. Make sure that the product is dry, if necessary dry the possible water residues with pressured medical air.
- Use the device and the autoclave appropriate bags of sterilisation in conformity with the NF EN868 standard. Immediately store the device after cleaning. Always store only a device per bag. According to the conditions of use, a double packing may be necessary.
- In order to avoid a retention of water, place the bag inside the autoclave so that the hollow parts are directed downwards.
- If the autoclave proposes several type of cycles, select a cycle for medical devices (at least 134°C, 2 bar, during 18 minutes.)

⚠ After sterilisation, check the absence of residual water inside and outside the packaging. Check the correct change of colour on the indicators.

- Keep the devices stored in sterilisation bag away from the light, the humidity and contamination of all kinds. Follow the recommendations of the manufacturer of packaging.
- The time of conservation of the device after sterilisation should not exceed 1 month.
- Label the devices by specifying the expiry date. Beyond the expiry date redo the cleaning and sterilisation cycles.

IX. REPAIR

In case of breakage, contact your authorised dealer or the ANTHOGRY After-Sales Service directly (+33 (0)4 50 58 50 53).

All the repairs must be carried out by manufacturer made parts and sub-parts. Repairs must only be done by one authorised dealer or by the After Sales Service of the manufacturer.

For any revision or repair, the medical device must be returned back, complete and **sterile**, with a certificate of sterility.

It must be joined with a document describing the technical problem while mentioning the name and address of the dental surgeon.

So that the requests for guarantees may be taken into consideration, please have the device returned with a copy of the invoice or of the delivery slip.

Replacement parts are available for 7 years beyond sales discontinuation.

X. ACCESSORIES

To be ordered from your authorised dealer.

Description	Reference
6 single use plastic filters	12006
12 single use plastic filters	12012
Adaptator Aspeo®/Aspiration saliva	12041

XI. WARRANTIES

This device is guaranteed parts and labour against all manufacturing defects for 12 months from the date of invoice. This guarantee does not apply to wear and tear parts. All changes or additions to the product without the express agreement of Anthogyr render this guarantee null and void. The guarantee becomes null and void if the technical instructions are not followed. Anthogyr cannot be held responsible for damage resulting from or which could result from normal wear, use, cleaning or incorrect maintenance, the non-observance of instructions for use or connection, scaling or corrosion, impurities in the water supply system or unusual chemical or electrical influences or non observance of the instructions, maintenance instructions and assembly of Anthogyr and other manufacturer's instructions. Delivery charges incurred when sending an instrument back to Anthogyr for repair will be paid by the client, even if the repair itself is covered by the guarantee. Postage and packing fees when returning the instrument to the client are covered by the guarantee. So that guarantee requests are taken into consideration, please attach a copy of the invoice or a copy of the delivery slip to the device.



2237 avenue André Lasquin
74700 SALLANCHES - FRANCE
Tél. +33 (0) 4 50 58 02 37
Fax +33 (0)4 50 93 78 60
N°SAV / Repairs : 33 (0)4 50 58 50 53
E-mail : sales@anthogyr.com
www.anthogyr.com

CE 0459