


axiom[®]

MULTI LEVEL[®]



Service Set
axiom[®]BL—TL

We would like to thank you very much for your trust and confidence as expressed in working with the Axiom® Multi Level® implant system.

Anthogyr offers various service solutions designed to solve problems encountered with the Axiom® Bone Level and Axiom® Tissue Level implant systems.

Success for you means success for us. Our marketing network and team of experts is completely available to you for any further information you may need.

Anthogyr



This instruction for use is not valid for OPKITRET, INKITEXPR kits.

These must be used according to the instructions for use available on the site ifu.anthogyr.com with search code: OPKITRET or INKITEXPR

INSTRUCTIONS AVAILABLE ONLINE

ifu.anthogyr.com

You can now find instructions for use (instructions and manuals) for **Anthogyr implants and prosthetics parts** in PDF format on our site ifu.anthogyr.com **using a PDF reader** (type Adobe Player).



HOW DOES THE SITE WORK?

This portal provides the latest instructions for using Anthogyr products.

To find the instructions for your device, please follow these steps:

1- Enter your product reference number, description or GTIN code (Global Trade Item Number) in the search field.

2- Press submit.

Your product's instructions will be available in PDF format, which you can consult online and/or print.

3- Select a language

Our instructions are available in several languages. To select the language you need, click the language choice menu.

This site is optimized for a 1024 x 768 px resolution screen to display instructions on PC or Mac with the following browser versions: Microsoft Internet Explorer 11 or higher, Safari 7.0 or higher (Mac only), Chrome 43 or higher, Firefox 38.0 and higher, and IOS and Android.

INFORMATION UPDATES:

Instructions for use are updated regularly and indicated by the "New" pictogram. Updated instructions can impact patient safety.

For this reason, we suggest you to avoid local back-ups and advise you to always check the Anthogyr portal.

To access archived instructions, click on "View old document versions".

You can also receive paper copies of instructions at no additional cost.

To receive paper copies, fill out the form available under the "Contact" tab or include a request with your next order.

Make sure to include your desired language.

The document will be delivered to you within 7 calendar days.

We are available if you have any comments or suggestions, via the "Contact" tab.

→ Warnings and recommendations

The instructions contained in this document describe the different phases of the surgical procedure and prosthetic restoration to be followed for the Axiom® Multi Level® implant system. A few general features specific to inserting implantable devices are recalled for information. This is not in any way an exhaustive document about implant and prosthetic practices to which the reader has any right of complaint.

TRAINING:

Axiom® Multi Level® components should only be used by practitioners who have been trained in implant practice and/or prosthetic techniques, and who are equipped for this type of procedure. Correct knowledge of surgical techniques and prosthetics is required to use this system. Specific training is offered and delivered at the Anthogyr company.

The surgical and prosthetic technique for the Axiom® Multi Level® system is performed exclusively in conjunction with the original components and instruments in accordance with the manufacturer’s recommendations. Anthogyr can take no responsibility in case of placement non-compliant with this manual and in case of use of implants or prosthetic parts or instruments foreign to the system.

Clinical evaluation of the patient and the choice of treatment solution are the sole responsibility of the practitioner. The implant diameter and length must be determined beforehand by the dental practitioner, depending on the clinical situation. Patients should also be informed of potential risks associated with implanting this type of device: oedema, bruising, haemorrhage, periodontal complications, transient or permanent nerve damage, local or systemic infections or inflammation, bone fractures, loosening or fracture of the implant, dehiscence, aesthetic problems, aspirating or swallowing the device, iatrogenic trauma, etc.



EQUIPMENT:

The practitioner using the system is responsible for the follow-up and maintenance procedures required to identify and treat any complications as early as possible and for ensuring the correct functioning and safety of the device. The references and the batch numbers of all components implanted, temporarily and/or definitively, must be recorded in the medical file of the patient.

Follow-up and maintenance are part of the knowledge of a practitioner trained in placing dental implants.

The practitioner is also responsible for defining the different settings for his/her equipment (instrument rotation speed, irrigation flow rate, etc), according to each clinical case, and for confirming that these are in good condition before each procedure.

Reusable instruments must be cleaned, decontaminated and sterilised before each surgery (even when first used) in accordance with current protocols in hospitals and clinics. The organisation of the operating room, preparation of operating staff and of the patient (premedication, anaesthesia, etc...) should follow current procedures and are the responsibility of the practitioner.

Anthogyr can under no circumstances be held responsible for any harm arising from defective handling or use. In order to avoid swallowing or inhaling small components, it is recommended that these are rendered secure by fixing them to the outside of the mouth with a suture thread. Whenever an instrument is changed, confirm that the contra-angle or key are correctly fixed by applying slight traction and ensure that each part is correctly fixed onto the transfer system outside the oral cavity.

CONSERVATION:

In producing our products, we have paid particular care and guarantee that a manufacturing control has been performed on all products made available for sale. In order to guarantee their integrity, it is recommended that they be stored in their original packaging at an ambient temperature of between 15 and 30°C, away from moisture and direct sunlight.

Protect packages from dust and do not store in the same premises as solvents and/or paints containing solvents or chemical substances. The device must be used before the expiration date indicated on the traceability label.

If the package (blister-closure / bag) is damaged or a defect is apparent when the product is opened, it is imperative that the device not be used and that the nature of the defect, part numbers and batch numbers of the components concerned are re-ported to the distributor or to Anthogyr. The technical specifications contained within these instructions are provided for indicative purposes only and cannot form the subject of any complaint.

The Axiom® Multi Level® implant system must not be used on animals.

Single-use devices must not be reused, or resterilised (risk of contamination and risk of alteration of functional surfaces).

The instructions for use here in may only be reproduced or disseminated with prior approval from the Anthogyr company. Anthogyr reserves the right to vary the technical feature of its products and/or to make changes or improvements to the Axiom® Multi Level® system without prior notice.

The Axiom® Multi Level® implant system is not compatible with other Anthogyr and competitors’ systems.

If uncertain, the user should contact the Anthogyr company before use.

This manual cancels and replaces all previous versions.

Explanations, symbols and diagrams			
	Device sterilized by Gamma radiation		Do not re-use, single-use device
	Manufacturing batch number of the device		Protect from light
	Commercial part number of the device		Do not use if packaging is damaged
	Manufacturing date of the device		Keep dry
	Expiration date of the device		Temperature restrictions
	Warning: observe instructions for use		Manufacturer
	Non-sterile device		Class I or Class IIa/IIb medical device complying with European Directive 93/42/CEE
	Serilise by autoclave outside of packaging		Tightening torque
	Do not use autoclave Sterilization		Global Trade item Number
	Complies with norms and standards in Russia		Federal (U.S.) law restricts this device to sale by or on the order of an authorised dentist
	Do not re-sterilise		See the precautions for use

SUMMARY TABLE

Step 1 Identify the implant		Step 2 Identify the problem encountered		Step 3 Solution to be applied	Step 4 Protocol to follow page number
oxiom® BL	The hexagonal impression is damaged			INMDS or INMDL	9
	The abutment is stuck in the implant but remains intact			INEXPS or INEXPL	12
	The M1.6 tapping of the implant is damaged			OPTAM16	10
	The M1.6 screw is broken in the implant			INKITSCWBL	14
	The inLink® lock is broken			INKITLOCK	18
	The abutment is broken in the implant	Single-component abutment	The abutment is broken at the M1.6 threading	INKITSCWBL	14
		Screw-retained abutment. Was the screw removed intact ?	YES	INKITABTBL	16
			NO	INKITABTBL and INKITSCWBL	16 + 14
oxiom® TL	The hexagonal impression is damaged			INMDS or INMDL	9
	The inLink® lock is broken			INKITLOCK	18
	The M2.8 tapping of the implant is damaged			ILTAM280	11
	The M1.6 tapping of the implant is damaged			OPTAM16	10
	The M1.6 screw is broken in the implant			INKITSCWTL	20
	The single-component abutment is broken in the implant at the M1.6			INKITSCWTL	20


The various KITS

Service solutions for Axiom® BL and TL implants are available depending on the analysis of the clinical case:


- Accurate knowledge of the problem: please refer to the **SUMMARY TABLE** to order the most suitable solution.
- Knowledge of the implant system only: order the Services SET(s) listed below.

⚠ All the KITS are supplied **non-sterile**, they must be cleaned and sterilised according to the cleaning and sterilisation manual available at ifu.anthogyr.com (search code: *INEXPS* for example).

1. KIT to be used in combination with an Axiom® BL implant




2. KIT to be used in combination with an Axiom® TL implant TL



3. KIT to be used in combination with a connection inLink®



INMDS / INMDL - S / L unscrewing mandrel



Indication:

Remove a prosthetic part / screw with a damaged hexagonal impression.

Composition:



S unscrewing mandrel
INMDS⁽¹⁾



L unscrewing mandrel
INMDL⁽¹⁾

⁽¹⁾For best performance and clinical results, we recommend not using these components for more than 10 times each.

Materials required in addition to the KIT:




Mandrel holder key
INCM


Protocol:

Step 1: Clean and sterilise the instruments.
Clean the hexagonal impression.

Step 2: Use the INMDS or INMDL unscrewing mandrel with the INCM manual wrench to hook the hexagonal impression to unscrew the screw: turn **counterclockwise** while applying pressure.



OPTAM16 - Tap M1.6 x 0.35



Indication:
Conform the damaged M1.6 tapping of an Axiom® BL or TL implant.

Composition:




Tap M1.6 x 0.35
OPTAM16⁽¹⁾

⁽¹⁾For best performance and clinical results, we recommend not using these components for more than 10 times each.

Protocol:


Step 1: Clean and sterilise the instruments.
Clean the hexagonal impression.

Step 2: While taking care to remain in the implant axis, tap gradually without forcing. After each half turn, unscrew the tap a quarter turn. Unscrew the tap regularly to clean it.



Step 3: Thoroughly clean the implant connector before inserting a new prosthetic part.

ILTAM280- Tap TL M2.8



Indication:
Conform the damaged M2.8 tapping of an Axiom® TL implant.

Must not be used with an inLink® abutment. If an M2.8 tapping on an inLink® abutment is damaged, replace the abutment.

Composition:




Tap TL M2.8
ILTAM280⁽¹⁾

⁽¹⁾For best performance and clinical results, we recommend not using these components for more than 10 times each.

Protocol:


Step 1: Clean and sterilise the instruments.
Clean the hexagonal impression.

Step 2: While taking care to remain in the implant axis, tap gradually without forcing. After each half turn, unscrew the tap a quarter turn. Unscrew the tap regularly to clean it.



Step 3: Thoroughly clean the implant connector before inserting a new prosthetic part.

INEXPS / INEXPL - Abutment Extractor S / L




Indication:

Remove an abutment stuck in an Axiom® BL implant (the abutment must not be broken).

If the abutment is broken, please refer to the **SUMMARY TABLE** to determine the service SET to be used.

Composition:




S abutment extractor
INEXPS⁽¹⁾



L abutment extractor
INEXPL



⁽¹⁾ Instrument supplied in the INMODOPP3 kit.




Protocol:

Step 1: Clean and sterilise the body and the rod.


Step 2: Assemble the rod in the body of the extractor leaving a gap of 2mm (approx.) between the head and the thumbwheel.




Step 3: Screw the instrument using the extractor body.



Step 4: Screw the thumbwheel until it presses against the body of the extractor to allow the removal of the abutment.



INKITSCWBL - KIT screw M1.6 BL broken



Indication:

- Remove a broken M1.6 screw (in the case of a screw-retained abutment) in an Axiom® BL implant.
- Remove a single-component broken abutment at the M1.6 threading in an Axiom® BL implant.

Composition:



Axiom® BL Ø1.80 service guide
OPGU180



Ø0.75 left helicoidal drill
OPFHG075⁽²⁾



mandrel holder key
INCM



unscrewing mandrel Ø0.75⁽¹⁾
INMD075

⁽¹⁾ For best performance and clinical results, we recommend not using these components for more than 10 times each.

⁽²⁾ This is a **single-use** component.

Equipment required in addition to the KIT:




manual key
OPCV060⁽³⁾-OPCV110-OPCV160⁽³⁾

⁽³⁾ These instruments are supplied in the INMODOPS3 kit.

Protocol:

Step 1: Clean and sterilise the instruments.
Clean the implant connector.
Lubricate the implant connector with water.

Step 2: Place the OPGU180 guide in the implant using an implant wrench (Ref. OPCV060, OPCV110, OPCV160) by exerting pressure to wedge it in the taper of the implant. Remove the implant wrench **while maintaining the guide in place**.




Step 3: Compress the OPGU180 guide.


⚠ Ensure that the guide is well anchored in the implant.

Lubricate the inner bore of the guide with water.

Step 4: Once the guide is in place, position the drill in the guide and then drill progressively (removing the shavings progressively) with the OPFHG075 drill **counterclockwise** at a speed of 3000 rpm until it stops. If the screw does not come with the drill, go to step 4.




Step 5: With the guide still in place, use the INMD075 unscrewing mandrel with the INCM wrench to hook the screw in the previously drilled Ø0.75 hole: turn **counterclockwise** while applying pressure.



Step 6: Thoroughly clean the implant connector before inserting a new prosthetic part.

INKITABTBL - Broken BL Pillar KIT




Indication:

Remove a broken screw-retained abutment in an Axiom® BL implant (abutment broken level with the neck of the implant).

If a single-component abutment is broken in an Axiom® BL implant, use the INKITSCWBL KIT (p. 14).


Composition:



retouch stop drill 1.6 mm
OPFBR16^[2]



Push rod + Extractor
INEXPRBL^[1]




stainless steel flat wrench 7 mm
INCP070

^[1] For best performance and clinical results, we recommend not using these components for more than 10 times each.


^[2] This is a **single-use** component.

Protocol:

Step 1: Clean and sterilise the instruments.
If the abutment is broken above the top of the implant, smooth out the surface of the abutment without damaging the implant.
If the abutment is broken below the top of the implant, do not adjust the abutment.
Clean the area with water and air.




Step 2: In the implant axis, under **irrigation**, drill with the OPFBR16 drill in a clockwise direction at a speed of 3000 rpm until it stops.




Step 3: Screw the extractor into the abutment while applying pressure. If required, finish tightening with the INCP070 wrench.

Ensure that the extractor is well anchored in the broken abutment.




Step 4: Screw the push rod as far as it will go in order to extract the abutment end attached to the extractor. If required, finish tightening with the INCP070 wrench.



Step 5: Thoroughly clean the implant connector before inserting a new prosthetic part.
If the M1.6 screw is broken in the implant, proceed with removing this screw using the INKITSCWBL KIT (p.14).

INKITLOCK - KIT broken lock



Indication:


Remove a broken inLink® lock when:

- The retaining ring and the head of the lock are in the prosthesis.
- The threaded part of the lock is in the Axiom® TL implant or the inLink® abutment.

Composition:



inLink® 2-in-1 wrench
IL-TOOL2



mandrel holder key
INCM



S unscrewing mandrel
INMDS⁽¹⁾



L unscrewing mandrel
INMDL⁽¹⁾

⁽¹⁾ For best performance and clinical results, we recommend not using these components for more than 10 times each.

Equipment required in addition to the KIT:



Spherical key
INBW100S⁽²⁾-INBW100L⁽²⁾


⁽²⁾ These instruments are supplied in the INMODOPS3 kit.

Protocol:

Step 1: Clean and sterilise the instruments.


Remove the threaded part of the lock:

Step 2: Use the INMDS or INMDL unscrewing mandrel with the INCM manual wrench to hook the threaded part of the lock: turn **counterclockwise** while applying pressure.



Remove the retaining ring and the head of the lock from the prosthesis:

Step 3: On the IL-TOOL2 wrench, remove the conical cap on the "IN" side.




Step 4: Insert the INBW100S ball wrench into the screw channel of the prosthesis until it engages in the head of the lock.




Step 5: Turn the prosthesis over so that the ball wrench head is pressed against the worktop to keep it engaged in the head of the lock.


Hold the prosthesis in one hand and engage the IL-TOOL2 wrench in the 2 spurs of the retaining ring in the other hand.



Step 6: Turn the IL-TOOL2 wrench clockwise to disengage the retaining ring and the lock head of the prosthesis.



INKITSCWTL - KIT screw M1.6 TL broken



Composition:



⁽¹⁾For best performance and clinical results, we recommend not using these components for more than 10 times each.

 ⁽²⁾This is a **single-use** component.

Equipment required in addition to the KIT:




⁽³⁾ These instruments are supplied in the INMODOPS3 kit.


Protocole:

Step 1: Clean and sterilise the instruments.
Clean the implant connector.
Lubricate the implant connector with water.


Step 2: Screw the TOPGU180 guide into the implant using a TIW100S implant wrench at a manual torque of <10N.cm until it stops. Remove the implant wrench.
Lubricate the inner bore of the guide with water.



Step 3: Once the guide is in place, position the drill in the guide and then drill progressively (removing the shavings progressively) with the OPFHG075 drill **counterclockwise** at a speed of 3000 rpm until it stops. If the screw does not come with the drill, go to step 4.






Step 4: With the guide still in place, use the INMD075 unscrewing mandrel with the INCM wrench to hook the screw in the previously drilled Ø0.75 hole: turn **counterclockwise** while applying pressure.








Step 5: Thoroughly clean the implant connector before inserting a new prosthetic part.

Component part numbers


AXIOM® BL KIT		REFERENCES
	BL M1.6 broken screw KIT Includes: - 1 Axiom® BL Ø1.80 guide - 1 unscrewing mandrel for Ø0.75 - 1 Ø0.75 left helicoidal drill - 1 mandrel holder key	INKITSCWBL
	BL broken abutment KIT Includes: - 1 retouch stop drill 1.6 mm - 1 push rod - 1 extractor - 1 stainless steel flat wrench 7 mm	INKITABTBL
	Axiom® BL service SET Includes: - the components of the INKITSCWBL KIT - 1 Axiom® BL Ø1.80 guide - 1 Ø0.75 left helicoidal drill - 1 unscrewing mandrel for Ø0.75 - 1 mandrel holder key - the components of the INKITABTBL KIT - 1 retouch stop drill 1.6 mm - 1 push rod - 1 stainless steel flat wrench 7 mm - 1 extractor - 1 S unscrewing mandrel - 1 L unscrewing mandrel - 1 S abutment extractor - 1 L abutment extractor - 1 M1.6 x 0.35 tap - 1 mandrel holder key	INKITSSBL

AXIOM® BL KIT		REFERENCES
	Broken lock KIT Includes: - 1 mandrel holder key - 1 S unscrewing mandrel - 1 L unscrewing mandrel - 1 inLink® 2-in-1 wrench	INKITLOCK

AXIOM® BL COMPONENTS		REFERENCES
	Unscrewing mandrel for Ø0.75 Medical Grade Stainless Steel Unscrewing mandrel for Ø0.75	INMD075
	Unscrewing mandrel Medical Grade Stainless Steel S Unscrewing mandrel L Unscrewing mandrel	INMDS INMDL
	Ø0.75 left helicoidal drill Medical Grade Stainless Steel Ø0.75 left helicoidal drill	OPFHG075
	M1.6 x 0.35 tap Medical Grade Stainless Steel M1.6 x 0.35 tap	OPTAM16
	Mandrel holder key Medical Grade Stainless Steel Mandrel holder key	INCM
	inLink® 2-in-1 wrench Medical Grade Stainless Steel inLink® 2-in-1 wrench	IL-TOOL2
	Stainless steel flat wrench 7 mm Medical Grade Stainless Steel Stainless steel flat wrench 7 mm	INCP070
	Axiom® BL Ø1.8 service guide Medical Grade Stainless Steel Axiom® BL Ø1.8 service guide	OPGU180
	Retouch stop drill 1.6 mm Medical Grade Stainless Steel Retouch stop drill 1.6 mm	OPFBR16
	Abutment extractor Medical Grade Stainless Steel S abutment extractor L abutment extractor	INEXPS INEXPL
	BL broken abutment extractor Medical Grade Stainless Steel BL broken abutment extractor	INEXPRBL

AXIOM® TL KIT		REFERENCES
	TL M1.6 broken screw KIT Includes: - 1 Axiom® TL Ø1.80 guide - 1 unscrewing mandrel for Ø0.75 - 1 Ø0.75 left helicoidal drill - 1 mandrel holder key	INKITSCWTL
	Axiom® TL service SET Includes: - the components of the INKITSCWBL KIT - 1 Axiom® TL Ø1.80 guide - 1 Ø0.75 left helicoidal drill - 1 unscrewing mandrel for Ø0.75 - 1 mandrel holder key - 1 stainless steel flat wrench 7 mm - 1 S unscrewing mandrel - 1 L unscrewing mandrel - 1 M1.6 x 0.35 tap - 1 TL M2.8 tap - 1 mandrel holder key	INKITSSTL

AXIOM® TL KIT		REFERENCES
	Broken lock KIT Includes: - 1 mandrel holder key - 1 S unscrewing mandrel - 1 L unscrewing mandrel - 1 inLink® 2-in-1 wrench	INKITLOCK

AXIOM® TL COMPONENTS		REFERENCES
	Unscrewing mandrel for Ø0.75 Medical Grade Stainless Steel Unscrewing mandrel for Ø0.75	INMD075
	Unscrewing mandrel Medical Grade Stainless Steel S Unscrewing mandrel L Unscrewing mandrel	INMDS INMDL
	Ø0.75 left helicoidal drill Medical Grade Stainless Steel Ø0.75 left helicoidal drill	OPFHG075
	M1.6 x 0.35 tap Medical Grade Stainless Steel M1.6 x 0.35 tap	OPTAM16
	TL M2.8 tap Medical Grade Stainless Steel TL M2.8 tap	ILTAM280
	Mandrel holder key Medical Grade Stainless Steel Mandrel holder key	INCM
	inLink® 2-in-1 wrench Medical Grade Stainless Steel inLink® 2-in-1 wrench	IL-TOOL2
	Stainless steel flat wrench 7 mm Medical Grade Stainless Steel Stainless steel flat wrench 7 mm	INCP070
	Axiom® TL Ø1.8 service guide Medical Grade Stainless Steel Axiom® TL Ø1.8 service guide	TOPGU180

NOTES

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

NOTES

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

 Anthogyr
2 237, Avenue André Lasquin
74700 Sallanches - France
Tél. +33 (0)4 50 58 02 37
Fax +33 (0)4 50 93 78 60
www.anthogyr.com

CE 0459