

CE 0459

Anthogyr
2 237, Av. André Lasquin
74700 Sallanches - France
Phone : +33(0)4 50 58 02 37
Fax : +33(0)4 50 93 78 60
www.anthogyr.com
Email : sales@anthogyr.com
Edition : 2016-03
REF. 0630PLA_NOT_EN



Medical device complying with
European Directive 93/42/EEC.



Attention: The accompanying
instructions must be consulted;
Danger.



See the precautions for use.



Do not reuse



Not sterile



Sterilisation by autoclave without
packaging, at the specified
temperature



Catalogue reference



Batch number



Do not use if the package is damaged



Store in a dry place with humidity
between 30% and 70%



Store away from light



Manufacturer



Manufacturing date



Federal (U.S.) law restricts this
device to sale by or on the order of an
authorised dentist.



Complies with norms and standards
in Russia

LOCATOR® Attachment System,

INCLUDES

Specific LOCATOR® Attachments, restora-
tive parts, dental laboratory components,
and instruments.

Important: This document contains the
most current Instructions for Use.
Please read and retain.

DESCRIPTION

Universal hinge, resilient attachment for
endosseous implants Axiom® REG or
Axiom® PX.

INDICATIONS

The LOCATOR® Implant Attachment Sys-
tem is designed for use with overdentures
or partial dentures, retained in whole or in
part, by Axiom® REG or Axiom® PX im-
plants in the mandible or maxilla.

CONTRAINDICATIONS

Not appropriate where a totally rigid con-
nection is required. Use of a single implant
with divergence of greater than 20 degrees
from vertical is not recommended.

CAUTION



Federal (USA) law restricts this de-
vice to sale by or on the order of a
licensed dentist.

WARNINGS, PRECAUTIONS & DIRECTIONS FOR USE

The LOCATOR® Attachment Systems have
not been evaluated for safety and compati-
bility in the MR environment. The LOCA-
TOR® Attachment Systems have not been
tested for heating or migration in the MR
environment.

SINGLE-USE DEVICES

The LOCATOR® Attachment Systems are
single-use devices. Single-use devices
must not be reused, not resterilized (risk of
contamination and risk of alteration of func-
tional surfaces)

If re-used, Anthogyr cannot guarantee the
functionality or the safety of the product.
Anthogyr accepts no responsibility for com-
ponents re-sterilized by the user.

STERILIZATION



All components and instruments are
supplied **NON-STERILE**. Titanium
abutments and instruments may be steri-
lized by Autoclave or Dry Heat sterilization
using the following parameters:

1. Autoclave sterilize using 121° C (250 ° F),
(15-20 psig at sea level), for 20 minutes
minimum.
2. Dry Heat sterilize using 170° C (338° F)
for 2 hours minimum.

The nylon males may be sterilized/

disinfected using a liquid chemical sterilant.
In order to ensure that the nylon males are
sterilized/disinfected (all microorganisms
including Clostridium sporogenes and Ba-
cillus subtilis spores are eliminated), the
nylon males must be soaked for a mini-
mum of 3 hours in the liquid sterilant at
room temperature.

PROSTHETIC PROCEDURES AND TECHNIQUE MANUAL

Using a calibrated torque wrench, tighten
the Locator Abutment to 25 Ncm



Warning: Use of higher torque values
than recommended above could
cause a fracture of the LOCATOR® Abut-
ment and/or implant.

The Technique Manuals for the LOCATOR®
Attachment System, are available from the
manufacturer or your distributor. These are
also available online at
www.zestanchors.com

The instructions are not sufficient by them-
selves for risk-free implantation of Axiom®
REG and Axiom® PX' implants. You abso-
lutely must follow the instructions of the
surgical manual "Axiom® REG and Axiom®
PX user guide". This document is delivered
with each surgical kit. It is available at
www.anthogyr.com or by simply requesting
it from Anthogyr at the contacts above.