

SAFE



IMPACTEUR AUTOMATIQUE

AUTOMATIC IMPACTOR

AUTOMATISCHER IMPAKTOR

IMPACTADOR AUTOMÁTICO

IMPATTATORE AUTOMATICO

IMPACTADOR AUTOMÁTICO

AUTOMATISCHE IMPACTOR

自动敲击器

آلة سحق بالصدم أوتوماتيكية

FR - NOTICE D'INSTRUCTION

EN - INSTRUCTIONS FOR USE

US - INSTRUCTIONS FOR USE

DE - GEBRAUCHSANWEISUNG

ES - INSTRUCCIONES DE USO

IT - ISTRUZIONI PER L'USO

PT - INSTRUÇÕES DE USO

NL - GEBRUIKSAANWIJZING

ZH - 使用说明

AR - تعليمات الاستخدام

TABLE OF CONTENTS

- I. DESCRIPTION OF SYMBOLS USED
- II. INDICATIONS FOR USE
- III. AREA OF APPLICATION
- IV. GENERAL SAFETY INFORMATION
- V. TECHNICAL SPECIFICATIONS
- VI. MD START UP PROCEDURE
- VII. HYGIENE AND MAINTENANCE
- VIII. REPAIRS
- IX. GUARANTEE
- X. ACCESSORIES - INSTRUMENTS
- XI. CONDITIONS OF STORAGE & TRANSPORT
- XII. DISPOSAL OF THE PRODUCT

I. DESCRIPTION OF SYMBOLS USED

	DANGER		Sterilisable up to 135°C
	Please wear gloves		Heat disinfectable
	Visual examination		Manufacturing date
	Please read this!		Manufacturer
	Batch number of the device		Medical Device (Impactor + Instruments)
	Device reference		Non-sterile device

II. INDICATIONS FOR USE

The **MD** (Impactor + Instruments) is intended for dental surgery. The **MD** impactor is designed to generate shocks with a controlled intensity and duration. The number of micro-impacts is regulated by the motor speed.

It is intended for use in the following areas:

→ Osteotomy:

Combined with the "osteotome" instruments, this **MD** impactor will facilitate the performance of the "tapping" type osteotome technique without the need for a mallet. It will enable sinus lift procedures to be performed using the crestal approach and osseous condensations in the jaw where bone mineral density is low.

→ Tooth extraction:

Combined with the "periosteal" instruments, this **MD** impactor will initiate tooth mobility by enlarging the desmodontal space in the interproximal, lingual and distal areas with the ultimate aim of atraumatic extraction. The use of periosteal helps preserve the alveolar bone.

→ The prosthesis:

Combined with the "end-piece impaction" instruments, this **MD** impactor will facilitate the placement of prosthetic components by means of impaction for the Axiom 2.8 range.

III. AREA OF APPLICATION

This **MD** (impactor + instruments) is intended exclusively for professional use in the field of dental surgery. Inappropriate or incorrect use could damage this device and result in hazardous risks to the user and to third parties.

In accordance with these provisions, the **MD** (impactor + instruments) may only be used by a user with experience in dental medicine, for the described application, and in line with applicable provisions concerning the prevention of accidents at work, and protection at work, as well as the indications in these operating instructions. Preparation and maintenance of this **MD** (impactor) may only be carried out by persons who have been trained in the field of infection control, in self-protection and in the protection of patients.

- In accordance with these provisions, the user must:
- Use working instruments solely without defects,
 - Observe the correct intended use
 - Protect themselves, as well as the patient or third parties, against any dangers
 - Prevent any contamination by the product

The following situations:

- Inappropriate use,
- Lack of maintenance,
- Use of accessories, instruments and spare parts that have not

been authorised by Anthogyr.

- Use of accessories and instruments from this **MD** with other devices.

- Modifications or additions not validated by Anthogyr to **MD** relieve Anthogyr of any obligation with regard to the guarantee or other claims.

This **MD** comply with Community Directive 93/42/EEC, amended by Directive 2007/47/EEC.

IV. GENERAL SAFETY INFORMATION

Before use, ensure that the **MD** (Impactor + Instruments) is not damaged and that there are no parts missing. Check the condition of the instruments used and handle them with due care and attention. The practitioner should test the **MD** Impactor without the patient to ensure that it is working properly.

Protective gloves should preferably be used whenever the tools are handled. Wear suitable protective clothing, particularly gloves, a mask and goggles.

After insertion of an instrument, check its correct retention with a slight axial movement.

Resterilisation of reusable **MD** must be performed by properly trained and protected personnel, in accordance with the applicable regulations.

In the case of visible malfunctions or damage, stop using the instrument immediately and inform your authorised distributor or the manufacturer.

In the case of any questions regarding the device, contact either your authorised distributor or the manufacturer.

All repairs must be carried out using parts and subassemblies certified by the manufacturer.


The use of a rubber dam in order to prevent tools falling into the mouth is strongly recommended.

Only use **MDs** (Instruments) that are compatible with the **MD** (Impactor).

Do not use instruments from this device on any other type or brand of **MD** or with any other device.

V. TECHNICAL SPECIFICATIONS

REFERENCE	69200ST/6920
Total length (mm)	142
Sleeve length (mm)	124
Maximum sleeve diameter (mm)	20
Rod diameter (mm)	7
Weight (g)	195
Connection tool	Standardised mandrel (ISO 1797-1) Mandrel with tri-lobe (proprietary Anthogyr connection) Locking through rotation of the ring
Connection motor	ISO 3964
Maximum torque transfer with the MD	Using the "tri-lobe" tools: 200 N.cm Using the ISO-standard tools: 80 N.cm
Energy	Mechanics [impacts]



Relationship between the input rotational speed (motor) and the output impact frequency

Ratio 1:1 rpm	Safe Impactor Impacts per minute
2 000	20
5 000	50
10 000	100
MAX 15 000	150

VI. MD START UP PROCEDURE

 The MD (Impactor + Instruments) are supplied **non-sterile** and **non-lubricated***.

Before using for the first time, the MD (Impactor + Instruments) must be **cleaned, decontaminated, lubricated* and sterilised** (see the section on "Hygiene and maintenance").

*: for the "Impactor SAFE3" MD only

 Resterilisation of reusable MD (Impactor + Instruments) must be performed by properly trained and protected personnel, in accordance with the applicable regulations. The resterilisation protocol must be adapted for a risk of infection.

 Before use, check that the MD (Impactor + Instruments) and/or the osteotomes are not damaged and that no parts are missing. Check the condition of the instruments used and handle them with caution and care. *(The list of accessories is included in paragraph IX)*


6.1. Motor connection:

 Protective gloves should preferably be used whenever the tools are handled.

Check the condition of the MD Impactor used and handle it with due care and attention.

 The MD (Impactor) may only be connected to the motor outside of the patient's mouth and when the motor has been stopped.

Position the MD (Impactor) on the motor coupling (in accordance with EN 23964 - ISO 3964) until it clicks into place. Keep the motor and the MD (Impactor) on the same axis for this purpose.



For certain motors, it is necessary to activate the side button to facilitate the connection.

Pull gently on the MD (Impactor) every time before use in order to verify that it is correctly engaged on the motor coupling.

Functional test

 Run the motor in the normal direction (the device will not operate if the motor is in "reverse" mode). Check that the MD (Impactor) properly generates regular micro-impacts.

NOTE: If you notice overheating, irregularities, vibrations or abnormal sounds during operation of the MD (Impactor), immediately contact your authorised technical service centre.

6.2. Disconnecting the motor:

 CAUTION: The MD (Impactor) may only be disconnected when the motor has been stopped, and outside of the patient's mouth. Withdraw the MD (Impactor) while maintaining the axis of the motor. For certain motors, it is necessary to activate the side button to enable the MD (Impactor) to be disconnected.

NOTE: To prevent oil from leaking into the motor, do not leave the MD (Impactor) connected to the motor if the MD (Impactor) is not used for prolonged periods.



6.3. Connecting / disconnecting the instrument:

 The instruments may only be connected to the MD Impactor when the motor is switched off and outside the patient's mouth.


 Protective gloves should preferably be used whenever the tools are handled.

 Check the condition of the instruments used and handle them with caution and care.


6.4. Locking/unlocking the MD (Instruments) in the MD (Impactor):

Positioning a tool:


- Align the arrow on the ferrule with the open padlock 




- If the MD (Instruments) has a trefoil symbol, align the symbol on the MD (Instruments) with the symbol at the end of the MD Impactor connection and insert the MD (Instruments).



- If the MD (Impactor) is standard, insert the tool by locating the drive bar.
- Turn the collet to align the arrow with the closed padlock and lock the MD (Impactor) on the MD (Impactor).




- Check that the MDs (Instruments) are firmly in place by applying gentle traction to them.



Removing a tool:

- Align the arrow on the ferrule with the open padlock



- Withdraw the tool by pulling on it

After insertion of an instrument, check its correct retention with a slight axial movement. It is essential that the motor is stopped and outside of the patient's mouth whenever the instrument is inserted or extracted.

The MD (Impactor + Instruments) does not require maintenance work. Nevertheless, should you notice a functional defect in the locking or gripping mechanism, this could mean that it is damaged. In this case, immediately return the device to the Anthogyr after-sales service department.

Make sure that the device does not have any areas of corrosion or cracks, and check that it is working properly.

MD lifecycle

If used in a proper manner, all MD parts have a lifecycle corresponding to 250 sterilisation cycles. However, these indications are not a warranty because wear may appear prematurely, depending on how the MD is maintained (cleaning and sterilisation).

VII. HYGIENE AND MAINTENANCE

Sterilisation of MD (Impactor + Instruments) must be performed by properly trained and protected personnel, in accordance with the applicable regulations. The sterilisation protocol must be adapted for the risk of infection.

In order to avoid any risk of infection and injury, it is essential to wear protective gloves.

Refer to the manufacturer's instructions for each product used. Observe the concentrations, lengths of exposure and the lifespan of the products. Do not combine products and observe the provisions regarding their disposal.

Only use products that are intended for the maintenance of medical and surgical equipment, compatible with stainless steel (no chlorine compounds), and compatible with aluminium (no soda).

Avoid the use of antiseptics that are intended for use on the skin and on mucous membranes.

Avoid the use of products containing aldehydes, alcohol or other products that are likely to bind proteins.

Preference should be given to detergents / disinfectants that are pH neutral or slightly alkaline.

The device must be cleaned and sterilised after every procedure.

The motor must be disconnected or stopped before handling.

Instruments must not be connected to the MD Impactor.

The MD Impactor is not removable and is not submersible.

The instruments are submersible.

7.1. Manual preparation for sterilisation:

Reprocessing of the MD (Impactor + Instruments) must be carried

out as soon as possible after treatment (within max. 2 hours).

By brushing:

- Brush under running water using a soft brush.
- Clean the entire MD (Impactor + Instruments) using a disinfectant wipe.

Or by spraying:

- Spray the disinfectant on to the MD (Impactor + Instruments) and wipe using a clean cloth.

Rinsing and drying

- Rinse properly¹
- Dry the entire MD (Impactor + Instruments).

7.2. Automatic preparation for sterilisation:

Only use a washer-disinfector appropriate for the reprocessing of this type of MD (Impactor + Instruments). Follow the instructions provided by the manufacturer of the unit.

- If the unit is equipped with nozzles for contra-angles: Fit the MD Impactor onto the nozzle of the device.

- If the unit is not equipped with nozzles for contra-angles: without disassembling the MD Impactor, position it with the tool connection component pointing upwards. Position the parts so that any retention of fluid is avoided. Immobilise each part.

- The thermal disinfection cycle must be at least 10 minutes at 93°C (203°F).

- Follow the instructions for use provided by the thermal disinfector manufacturer, and in particular the recommendations regarding the products that should or should not be used with the unit.

- Observe the drying cycle (Caution: Do not exceed 140°C (284°F) during the drying cycle)

- Check that there is no residue and that all parts of the MD (Impactor + Instruments) are perfectly dry at the end of the cycle.

The MD Impactor is not submersible.

7.3. Lubrication:

- Wear a protective face mask.

- Disconnect the MD (Impactor) from the motor before handling.
- Remove the connected instrument in the MD (Impactor)
- Check that the nozzle of the spray is the correct one for use with the MD (Impactor).
- Insert the aerosol spray at the back of the MD (Impactor).
- Cover the head (connecting part) of the MD (Impactor) with a wipe.
- Position the head (connecting part) of the MD (Impactor) downward.
- Spray several times, keeping the MD (Impactor) closed.
- Wipe off surplus oil with a soft cloth or a wipe.

- Lubricate after each decontamination and before each sterilisation, and repeat in the event of prolonged use.

- Keep away from heat and ignition sources. Smoking is strictly prohibited!

7.4. Functional test:

This test must be performed every time prior to sterilisation. Connect the MD (Impactor) to the motor and point the head downwards. Run at low speed (5000 motor rpm) for 30 seconds in order to clear the excess oil, then at full speed (10,000 motor rpm) for 30 seconds. Wipe the MD (Impactor) with a disinfectant wipe if the oil has dripped.

Note: If you notice overheating, irregularities, vibrations or abnormal sounds during operation of the MD (Impactor), immediately contact your authorised technical service centre.

7.5. Sterilisation:

Resterilisation of reusable MDs must be performed by properly trained and protected personnel, in accordance with the applicable regulations. The resterilisation protocol must be adapted for a risk of infection.

The MD (Impactor + Instruments) and the accessories must be sterilised before they are used for the first time, and every time after they are used (after cleaning and disinfection). Do not sterilise instruments unless they have been pre-disinfected, cleaned, decontaminated, lubricated (only the MD Impactor) and tested.

No disassembly of the MD is envisaged prior to performing sterilisation.

¹ Temperature <30°C. demineralised water is recommended instead of running water, for example if the latter contains too much chlorine (see the FDI89-135 standard)

We strongly advise the use of a class B autoclave.
All other methods of sterilisation must be avoided.
Read the instructions for use provided by the autoclave manufacturer.

Maintain the space between the bags and do not overload the autoclave.

Only sterilise instruments that have been predisinfected, cleaned, lubricated and tested.

Remove the instruments from the MD Impactor before sterilisation (See Section 5.5).

Make sure that the MD (Impactor + Instruments) do not have any areas of corrosion or cracks, and check that it is working properly. Ensure that the MD (Impactor + Instruments) are dry; if necessary, dry any residual water using medical quality pressurised air.

Use the sterilisation bags that are designed for the MD (Impactor + Instruments) and for the autoclave, or use a sterilisation cassette.

Each bag should always only contain one MD or sterilisation tray.

 In order to avoid any retention of water, place the bag in the autoclave in such a way that any hollow parts are face down.

If the pre-vacuum autoclave offers a number of different cycles, choose a cycle for MD **at least 135°C at 2.13 bar [275°F at 30.88 psi] for 3 minutes with a drying time of 16 minutes.**

After every sterilisation cycle, check that there is no water remaining on the inside and outside of the packaging.

Check that the change in colour of the flow indicator is correct.

7.6. Storage:

 Keep the MD (Impactor + Instruments) in the sterilisation bags away from light, moisture and contamination of any kind. Follow the recommendations of the packaging manufacturer.

The duration for which the MD (Impactor + Instruments) may be kept after sterilisation may not exceed 1 month.

Label the MD (Impactor + Instruments), specifying the expiration date. After the expiration date, repeat the cleaning and sterilisation cycle.

VIII. REPAIRS

In the case of device failure, please contact your authorised distributor, or contact our after-sales service department directly. Repairs may only be carried out by an authorised repair service or by the Anthogyr after-sales service department, and only using original Anthogyr replacement parts.

For any servicing or repair work, the MD must be returned complete and **sterile** with proof of sterility. It must be accompanied by a document that describes the problem in question and also includes the complete contact details of the practitioner who used the device.

Exchange of replacement parts is possible for 7 years after sales discontinuation.

CONTACT DETAILS: AFTER-SALES SERVICE

AFTER-SALES SERVICE DEPARTMENT

Anthogyr

2237 Avenue André Lasquin - 74700 Sallanches - FRANCE

Direct line: +33 (0)4 50 58 50 53

Mail: sav@anthogyr.com

IX. GUARANTEE

These MD (Impactor + Instruments) are guaranteed for parts and labour against any defect in manufacture for a period of **12 months** from the date of invoicing. This guarantee does not apply to wear parts and does not cover transport costs. For claims to be considered under the guarantee, please enclose a copy of the invoice or the delivery note with the MD (Impactor + Instruments). Any modification of or addition to the product without the express consent of Anthogyr will nullify this guarantee.

The guarantee shall become null and void if the technical instructions supplied with all of our equipment are not observed.

Anthogyr cannot be held responsible for damage and ensuing damage resulting or possibly resulting from normal wear and tear, incorrect use, cleaning or maintenance, non-compliance with the provisions regarding use or connection, scaling or corrosion, impurities in the water supply system or abnormal chemical or electrical effects, or non-compliance with the instructions for use, maintenance or installation provided by Anthogyr, and other manufacturer instructions.

X. ACCESSORIES - INSTRUMENTS

Instruments		
References	Designation	
OSTSCC34	Osteotome Ø3.4 concave straight	OSTEO SAFE Solution
OSTSCC40	Osteotome Ø4.0 concave straight	
OSTSCC46	Osteotome Ø4.6 concave straight	
OSTSCC52	Osteotome Ø5.2 concave straight	
OSTSCX34	Osteotome Ø5.2 convex straight	
OSTSCX40	Osteotome Ø4.0 convex straight	
OSTSCX46	Osteotome Ø4.0 convex straight	
OSTSCX52	Osteotome Ø5.4 convex straight	
OSTECC34	Osteotome Ø3.4 concave bayonet	
OSTECC40	Osteotome Ø4.0 concave bayonet	
OSTECC46	Osteotome Ø4.6 concave bayonet	SAFE LOOK Solution
OSTECC52	Osteotome Ø5.2 concave bayonet	
OSTECX34	Osteotome Ø3.4 convex bayonet	
OSTECX40	Osteotome Ø4.0 convex bayonet	
OSTECX46	Osteotome Ø4.6 convex bayonet	
OSTECX52	Osteotome Ø5.2 convex bayonet	
OPIPI100	Tip for abutment 0°/7°	
OPIPI200	Tip for abutment 15°/23°	
OPIPI400	Tip for extra-orally cemented prosthesis	
PERIO01	Cx right Fine Periotome	EXO SAFE Solution
PERIO02	Cx right Wide Periotome	
PERIO03	Cx Ang. Vertical Periotome	
PERIO04	Cx Ang. Horizontal Periotome	
PERIO05	Cc right Fine Periotome	
PERIO06	Cc right Medium Periotome	

Accessories

References	Designation	
1930X	Lubricating spray	All solutions
1932X	ISO type E connection tip	

XI. CONDITIONS OF STORAGE & TRANSPORT



XII. DISPOSAL OF THE PRODUCT

The MD (Impactor + Instruments) must be sterilised before disposal in order to prevent the risk of third-party contamination. On the basis of current knowledge, the product does not contain any substances that are harmful to the environment. Observe national legislation, standards and provisions with regard to disposal.

2237 avenue André Lasquin
74700 SALLANCHES - FRANCE
Tél. +33 (0) 4 50 58 02 37
Fax +33 (0) 4 50 93 78 60
N°SAV / Repairs :
+33 (0) 4 50 58 50 53
E-mail : sales@anthogyr.com
www.anthogyr.com